

UNITY-xr

The most versatile and powerful thermal desorption unit












Introducing the UNITY-xr thermal desorber for GC and GC-MS – the most versatile instrument for the analysis of trace-level volatile and semi-volatile organic compounds (VOCs and SVOCs) in air and materials.

As the world leader in thermal desorption (TD) since 1997, Markes International has pioneered many breakthroughs in analytical instrumentation. We now present the UNITY-xr for method-compliant TD analysis, across a wide range of sample types and applications.

Together with the other members of the new 'xr' series, the UNITY-xr offers




Outstanding analytical capability and reliability

Leading-edge thermal desorption technology


Building on the field-proven strengths of its UNITY 1 and UNITY 2 predecessors, every UNITY-xr offers:

- Analysis of ultra-volatiles to semi-volatiles: Optimised temperature capabilities, an upgraded flow path, and a new optional module for water management allow quantitative recovery of C₂ to C₄₄, including reactive and thermally labile species, from percent to sub-ppt concentrations.
- Sample splitting and re-analysis: Quantitative re-collection of both the inlet (tube desorption) split and the outlet (trap desorption) split is offered as standard on UNITY-xr (and all other TD systems from Markes), for repeat analysis and powerful method validation.
- Outstanding reliability: High-quality components reduce the need for routine maintenance, so maximising uptime.

Method-compliant air monitoring



Analysis of trace-level VOCs in ambient air, in accordance with US EPA TO-17 or Chinese EPA Method HJ 644, is easy with the UNITY-xr. Splitless analysis ensures maximum sensitivity for small sample volumes.



Unbeatable application versatility

Fully method-compliant thermal desorption analysis of sorbent tubes

The enhanced features of the UNITY-xr, coupled with a suite of innovative sampling accessories, make it suitable for a wide range of applications. Across many of these areas, our involvement with technical committees and legislative agencies means that we are uniquely well-placed to advise on method compliance.

Environmental monitoring



The UNITY-xr complies with:

- US EPA Method TO-17 (ambient air)
- US EPA Method 325 (fenceline)
- Chinese Method HJ 644 (ambient air)
- New European SVOC protocols
- Method CEN/TS 13649 (stack emissions)
- Chinese Method HJ 734 (source emissions)
- and more...

Indoor and in-vehicle air



The UNITY-xr complies with:

- ISO 16000 series (indoor air)
- ASTM D6196 (indoor air)
- ISO 12219 series (automotive test)
- VDA 278 (automotive test)
- Multiple OEM standards
- and more...


Consumer environmental health




The UNITY-xr complies with:

- European Method EN TS 16516 (construction products – mandatory)
- ASTM standards (spray polyurethane foam)
- New CEN standards (combustible air fresheners)
- and more...

Perfect for challenging applications





Defence and homeland security



Fragrance and odour profiling



Food and drink



Forensic



Biological profiling



Using Markes' dedicated PAH sorbent tube and a sample volume far lower than required by traditional monitoring protocols, the UNITY-xr allowed several ppt-level polycyclic aromatic hydrocarbons (PAHs) to be detected in this complex urban air sample.

Modular systems for ultimate sampling flexibility

The 'xr' series - For tubes, canisters, bags and on-line monitoring

The modular nature of Markes' systems means that it is easy to extend the analytical capability of any UNITY-xr, for even greater versatility.

CIA Advantage

Automated canister/bag autosampler

Added capability for:

- Up to 27 canisters/bags.
- Canister/bag monitoring of C₂-C₁₄.
- Manual re-collection of outlet split from canister/ bag samples.

Capability for:

Desorption of single tubes. Single-tube thermal desorber

UNITY-xr

- Monitoring of C₂-C₄₄ and reactive species.
- Manual re-collection of inlet and outlet split flows.

Air Server-xr

On-line VOC monitoring system

Added capability for:

- Up to eight channels.
- On-line monitoring of C_2 - C_{14} .
- Manual re-collection of outlet split from on-line samples.




Added capability for:

- 100-tube automation.
- Automated re-collection of outlet split flows.
- Automated re-collection of canister and on-line samples.
- Fully automated sequencing for tubes and canisters/bags.



Multi-tube autosampler



- Automated re-collection of inlet and outlet split flows.
- Automated re-collection of split flows from all 100 tubes onto 100 clean tubes.





High-throughput autosampler for UNITY-xr for up to 100 tubes. It offers:

- Automated re-collection options.
- Enhanced traceability of samples using barcodes and RFID TubeTAG.





- Stringent sample sealing throughout the desorption process using Markes' patented DiffLok caps.
- Additional confidence in results: Optional internal standard/dry-purge (ISDP)





Air Server-xr

10

15

Versatile on-line VOC
monitoring system which
can be used for on-line
analysis of ozone precursors,
terpenes, ultra-volatile
compounds, sulfur species,
greenhouse gases, in
addition to typical 'air toxics'.
It also complies with various
standard protocols, most
notably the US EPA PAMS scheme





25

30

35

20

Low limits of detection: Efficient operation of the UNITY-Air Server-xr provides LODs well below the 0.5 ppb required for analysis in accordance with the US EPA PAMS scheme, for ozone precursors ranging in volatility from acetylene to trimethylbenzene.


CIA Advantage

Cryogen-free automated canister autosampler and pre-concentrator.

Systems containing CIA Advantage-xr are versatile, offering gas-loop and MFC-controlled sampling options,



which combined with sample splitting, allow compatibility with the widest possible sample concentration range.



Method-compliant 'air toxics' analysis: This analysis of 1 L of a 1 ppb 'air toxics' mix shows how canister-sampled volatile air pollutants can be analysed in accordance with US EPA Method TO-15 or Chinese EPA Method HJ 759.

Markes Instrument Control

Easy-to-use software for the new 'xr' series




The new software used to control the UNITY-xr and the other members of the 'xr' series offers the following features for enhanced laboratory productivity:

- Automated, unattended sequencing of tube and on-line samples.
- **Editing of active sequences,** for greater flexibility and ease of use.
- Rapid set-up of TD methods using pre-programmed parameters for standard methods including VDA 278, US EPA TO-17 and PAH analysis.
- Pre-loading of an internal standard on a tube or trap, for enhanced quantitation.
- System self-checking, for improved diagnostics.

Unmatched product range

A comprehensive range of sorbent tubes and sampling accessories for every TD application



Markes International – The TD experts

World-leading instruments and unmatched expertise in VOC and SVOC monitoring

Markes International has for 20 years been at the forefront of innovation for enhancing the measurement of trace-level VOCs and SVOCs by thermal desorption-gas chromatography. Our suite of instruments for thermal desorption sets the benchmark for quality and reliability:

TD100-xr[™]High-throughput
100-tube automated thermal desorber.

UNITY-Air Server-xr[™] Versatile on-line VOC monitoring system. **ULTRA-xr[™]**High-throughput 100-tube autosampler for UNITY-xr.

CIA Advantage™ Cryogen-free automated canister autosampler and pre-concentrator. TC-20[™] & TC-20 TAG[™]
Cost-effective systems for off-line multi-tube conditioning and dry-purging.

TT24-7™
Twin-trap instrument for near-real-time on-line monitoring.

Micro-Chamber/Thermal Extractor™ Unique sampling device for emissions of VOCs and SVOCs from products and materials.

Markes International

UK: Gwaun Elai Medi-Science Campus, Llantrisant, RCT, CF72 8XL T: +44 (0)1443 230935

US: 11126-D Kenwood Road, Cincinnati, Ohio 45242

2355 Gold Meadow Way, Gold River, Sacramento, California 95670

T: 866-483-5684 (toll-free)

Germany: Schleussnerstrasse 42, D-63263 Neu-Isenburg, Frankfurt **T:** +49 (0)63

T: +49 (0)6102 8825569

