






# Sorbents – A quick reference guide

Sorbent	Volatility range		Suitable analytes	Max. temp (°C) [1]	Desorb temp (°C) [2]	
	Hydrocarbon range	Boiling point (°C)				
<b>Porous polymer</b> 	<b>Tenax TA</b>	C <sub>6</sub> ~C <sub>30</sub>	100 ~ 450	Aromatics, apolars, polar compounds bp >150°C, and semi-volatiles.	350	300
	<b>Tenax GR</b>	C <sub>6</sub> ~C <sub>30</sub>	100 ~ 450	Aromatics, apolars, polar compounds >150°C, and semi-volatiles.	350	300
	<b>HayeSep D</b>	C <sub>5</sub> ~ C <sub>12</sub>	50 ~ 200	Low molecular weight compounds, acetylene, halogens, and sulfur groups. GB/GE derivative of VX (CWA).	290	280
	<b>PoraPak N</b>	C <sub>5</sub> ~ C <sub>10</sub>	50 ~ 200	Polar VOCs, acrylonitrile, acetonitrile, propionitrile, pyridine, volatile alcohols, ethanol, methyl ethyl ketone.	190	165
	<b>PoraPak Q</b>	C <sub>5</sub> ~ C <sub>12</sub>	50 ~ 200	Oxygenated compounds.	250	190
<b>Graphitised carbon black</b> 	<b>Carbograph 2TD</b> <b>Carbopack C</b> <b>Carbotrap C</b>	C <sub>8</sub> ~ C <sub>20</sub>	130 ~ 340	Alkyl benzenes and large aliphatics. Heavy organics: PCBs, PNAs.	400	360
	<b>Carbograph 1TD</b> <b>Carbopack B</b> <b>Carbotrap B</b>	C <sub>5</sub> ~ C <sub>14</sub>	50 ~ 250	A wide range from medium to high volatility: Ketones, alcohols, and aldehydes (but not formaldehyde). Non-polars within volatility range. Perfluorocarbon tracer gases. BTX.	400	360
	<b>Carbograph 5TD</b> <b>Carbopack X</b>	C <sub>3</sub> ~ C <sub>8</sub>	50 ~ 150	Especially good for 1,3-butadiene and light hydrocarbons.	400	360
	<b>Carbopack Y</b> <b>Carbotrap Y</b>	C <sub>12</sub> ~ C <sub>20</sub>	50 ~ 340	Less volatile hydrocarbons. Pesticides e.g. alachlor atrazine, isoprene, and formothyon.	400	360
<b>Carbonised molecular sieve</b> 	<b>SulfiCarb</b>	C <sub>3</sub> ~ C <sub>8</sub>	-30 ~ 150	Thiols. WVOCs (vinyl chloride, CS <sub>2</sub> , methanol, ethanol, and acetone). Used for sterically large WVOCs (SF <sub>6</sub> ).	400	360
	<b>Carbosieve S-III</b>	C <sub>2</sub> ~ C <sub>5</sub>	-90 ~ 80	Permanent gases, and ethene to n-C <sub>5</sub> , e.g. chloromethane. Also ethylene from small volumes.	400	360
	<b>Carboxen 1000</b>	C <sub>2</sub> ~ C <sub>5</sub>	-60 ~ 80	Permanent gases, and ultra-volatile hydrocarbons, e.g. vinyl chloride.	400	360
	<b>Carboxen 1003</b>	C <sub>2</sub> ~ C <sub>5</sub>	-60 ~ 80	Permanent gases, and ultra-volatile hydrocarbons, e.g. ethane.	400	360
<b>Zeolite molecular sieve</b> 	<b>Molecular Sieve 5Å</b>	C <sub>2</sub> ~ C <sub>4</sub>	-90 ~ 80	Permanent gases, and nitrous oxide	400	300 (165 for N <sub>2</sub> O)
<b>Other</b> 	<b>Silica Gel</b>	N/A	N/A	Low-boiling polar compounds, especially useful for separating chlorinated or sulfur compounds from matrices with hydrocarbon interferences.	200	180
	<b>Quartz wool</b>	N/A	N/A	SVOCs.	400	360

[1] Maximum temperature should not be exceeded, to prevent irreversible damage to sorbent.

When using multi-beds of sorbents, maximum temperature is the lowest of all sorbents within the tube.

[2] Recommended desorption temperatures are a guide and may require optimisation depending on target analytes.

Lower desorption temperatures are recommended where possible to prolong sorbent lifetime. Sorbent conditioning should be carried out at 10–20°C above the desorption temperature.