


# Robotic Sample Preparation for Analytical Laboratories

Accurate, Affordable, and Adaptable Workflow Automation



## ePrep Sample Preparation Workstation

**ePREP**<sup>®</sup>  
automation for every laboratory

# ePrep® Sample Preparation Workstation

Designed to automate analytical syringe operation for ultimate precision and accuracy, ePrep helps you achieve:



## GREATER EFFICIENCY AND PRODUCTIVITY

Take 5-15 minutes to program your workflow and ePrep will do the rest for you



## SIMPLE METHOD OPTIMIZATION AND VALIDATION

ePrep's highly adaptable software allows you to test different variables to find the best performing method



## IMPROVED RESULTS

Eliminate human error and variability in sample preparation skill to easily validate your methodology



## GREENER PROCESSES

Automate with syringes allows you to accurately use smaller volumes of reagents and solvents resulting in less wastage

## ePrep Capabilities

---

### Precise Control of Small Volume Dispensing

ePrep boasts unparalleled accuracy and precision for liquid handling, capable of dispensing small volumes down to 0.1 µL

### “High Pressure” Applications

Suitable for applications requiring pressurised, positive displacement such as microSPE and filtering

### Suitable for Volatile Organics

With analytical syringe, sealed vials can be used to contain volatile organics to prevent sample loss and potential hazard. Important in applications such as VOC and toxic samples

### Versatile Robotic Syringe Change


Syringe volume can be changed during a workflow to reduce dispense error and improve efficiency (reduced washing)

### Thorough Tool Wash

ePrep washes syringes between tasks, effectively eliminating detectable carryover. Single use syringes can also be used



## Designed for Analytical Chemistry



### The Dilemma of Sample Preparation in Analytical Laboratories

Typically, an analytical laboratory has different sample preparation workflows [type (A) – (F) with varying sample batch sizes (5 – 185 samples) analysed on different instruments.

ePrep helps analytical laboratories by providing a flexible automation solution capable of switching rapidly between different methods even if for a small number of samples.

## Highly Flexible and Adaptable

- Rapid programming and modular configuration allow the ePrep to be set up and operational within minutes
- Independent operation allows a single ePrep to feed multiple analytical instruments uninterrupted, maximising efficiency

## Easy to Use

- ePrep's proprietary touch screen software is designed to be used independent of user skill. Workflows are created and executed in just minutes using 'drag-and-drop' style operation
- Features workflow validation checking, error feedback, and tool identification

## Complete Workflow

Using analytical syringes rather than pipettes, ePrep can precisely control micro volumes and flowrates at high pressure allowing separation, filtering, and membrane tasks to be used in a workflow sequence


## Target Applications

- Ideal for preparing samples of small-to-medium batch sizes, typically found in chromatography laboratories
- **Application Areas:** Especially suited for critical environmental, clinical, biotechnology, and food sample preparation using small volumes, sealed vials, and pressure applications with analytical syringes
- **Calibration Standard:** Excellent for preparation of samples with standard additions according to their associated reference calibration standards

## Validated Results


ePrep allows samples to be prepared the same way every time eliminating errors, saving time, and avoiding repeat:

- With typical volume accuracy and reproducibility of  $\leq 0.3\%$ , the improved data quality produced using the ePrep eliminates the need for duplicate and triplicate samples
- Predefined wash and purge routines minimise carry over eliminating errors



Calibration curve of Tri-nitrobenzene from an ePrep 5µL-80µL serial dispense.  $R^2 = 0.9997$

# Simple Setup




# Key Features

## COMPLETE AUTOMATION



Dilution, sample aliquot, standard addition, filtering, microSPE...and more. A range of add-on accessories are available for complete workflow automation

## FLEXIBLE CONFIGURATION




Rapid configuration of vial racks, accessories and elements to perform a range of sample preparation tasks

## WORKFLOW AND CONFIGURATION DEPOSITORY



Sample preparation methods and calibration standards can be stored for repeated use and sharing with other laboratories

## ROTATING FOOT




Enables syringe interchange, filter, and SPE cartridge handling

## ANALYTICAL SYRINGES




Analytical syringes enable high precision and accuracy dispensing down to 2µL

## ROBOTIC SYRINGE CHANGE




XCHANGE® technology enables automated and timely exchange of the necessary syringes and tools in a Workflow

## CLEANING STATION



Syringe wash station for zero carryover, clean sample processing, and waste disposal

## AUTOMATED MICROSPE & FILTERING



µSPEed cartridges enable high Resolution SPE for clean, low volume, no blowdown prep. Standard Disk Filters automate Sample filtering

## Tool Options

---



**ePrep Needle Syringes**

100µL, 1mL & 10mL  
for liquid dispensing



**ePrep eZy-Connect™ Syringes**

100µL, 1mL & 10mL  
for µSPEed cartridges



**ePrep Vial Gripper**

for transfer of vials and  
barcode read

## Popular Accessory Options

---




µSPEed Cartridge Rack




50mL Reagent Jar Kit



Bulk Solvent Manifold



Vial Vortex Mixer




Barcode Reader Module

## Rack/Vial Adapter Plates

---

To allow the processing of samples from sample vials directly into autosampler racks, ePrep has a range of deck Adapter Plates for most common sample and instrument rack designs. These include **BelArt, Supelco, Wheaton, 96 Well Plates, Agilent, Shimadzu, Thermo, Waters and CTC-PAL**. Each Adapter Plate has a corresponding software profile for simple programming of vial identification and position within an ePrep Workflow



# Technical Specifications

---

## Model

GL950X, P/N 01-01000-01, ePrep Sample Preparation Workstation

## LIQUID HANDLING AND SAMPLE PREPARATION

### Liquid Handling Precision

Liquid: 0.3% RSD dispensed at 5% syringe volume position.  
Mechanical:  $\pm 50\mu\text{m}$  plunger displacement precision.

### Liquid Handling Accuracy

$< \pm 1.0\%$  of total syringe volume with an uncalibrated syringe, Syringes can be calibrated for greater accuracy.

### Carry-over

For multi-use liquid handling tools an integrated wash station is used to control and eliminate carry over. The carry-over actually achieved depends on the syringe type, sample characteristics, and user set parameters in the software.

### Volume Range

Minimum volume: 5 $\mu\text{L}$  (precision of lower volumes is based on operating conditions).  
Maximum volume: 10mL per dispense with unlimited volume achievable using multiple aspirate/dispense cycles.

### Flow Rates

100 $\mu\text{L}$  syringe: 5nL/min - 18.0 mL/min  
1mL syringe: 50nL/min - 33.0 mL/min  
10mL syringe: 500nL/min - 36.0 mL/min

## OPERATIONAL INFORMATION

### Sample Capacity

Samples/Reagents can be contained in a wide variety of vial sizes and racks depending on deck configuration and application. Typical dispense volumes are between 10 $\mu\text{L}$  to 10mL.  
Total 880mm x 330mm deck space.  
Supports a wide range of vials from 96 well micro titer plates through typical 1.5mL vials to 60mL bottles.

### Laboratory

Sealed vial operation means a cover is not required. It is suggested well ventilated and climate-controlled environment is used. If using toxic materials, take appropriate safety precautions.

### Safety

Rapid emergency stop is activated when the ePrep system detects the presence of personnel and obstructions in the work area.

### Tablet Controller and Software

Microsoft® Touch Screen Surface Pro® (Supplied) including ePrep's Axis rapid Workflow development software (English only) and lifetime updates.

### Connectivity Tablet-to-Instrument

USB or Bluetooth

### Data Security

Multi-level rights controlled by user login.  
Software logs Date/Time, Output ID, Workflow ID, Processes. Log report formats in TXT and PDF.  
21CFR Part 11 compliance.

### Wash Station

Deck mounted syringe Wash Station. Pump driven with active waste removal. External solvent and waste reservoirs (reservoirs not supplied).

### Asset Identification

RFID identification on syringes, tools and racks. Optional 1D and 2D bar code identification on vials.

### Available Accessories (Optional)

Shaker, Gripper, Barcode Reader, Direct Detector Interface, Multiport Reagent Manifold,  $\mu\text{SPEed}$  Cartridge Rack

## SYSTEM INFORMATION

### Instrument Dimensions (L x W x H)

1370mm x 694mm x 743mm

### Weight

65kg

### Voltage

24V, 221W, 9.2A DC Power Pack; 100-240 Volt AC input; Requires country specific IEC C13 power cable (not supplied). Battery backup for continual operation even with power interruption.

### ePrep Connection Ports

3 x Serial Comm Ports, 4 x Digital Output, 4 x Digital Input, 4 x Relay

### Operating Temperature

10-35°C, 0-80% relative humidity

### Operating Sound Level

Typical 80dB

### Compliance

CE, FCC, IC, RCM, RoHS, Safety EN 61010, EMC 61326

### Warranty


12 months

*Specifications are subject to change without notice*

# Ordering Information

---

## ePrep Workstation

| Image   | Part No            | Description   |
|---|--------------------|---|
|  | <b>01-01000-01</b> | ePrep Sample Preparation WorkStation<br>Supplied with: Wash Station, Power Supply, Surface Pro Control Tablet, Cabling and Starter Kit including Adapter Plates, Vial Racks, Vials, ePrep Syringes and Spares |





## Adapter Plates

See Website for Sample and Autosampler Rack Adapter Plates

## Syringes and Tools

See Website for Syringes and Tools

## Accessories

| Image   | Part No            | Description  |
|---|--------------------|--|
| <u>Vortex Mixer</u>   |                    |  |
|  | <b>01-04100-01</b> | Vortex Mixer kit including power supply<br>Country-specific IEC C13 power cable is not included                          |
| <u>μSPEed Cartridge Rack</u>  |                    |  |
|  | <b>01-04160-01</b> | Rack for μSPEed Cartridges (96 cartridges)<br>Supplied with 500 μL ePrep Syringe<br><br>Note: Cartridges sold separately |
| <u>Barcode Reader Module</u>  |                    |  |
|  | <b>01-04010-01</b> | Barcode Reader Module<br>Supplied with gripper and installation tools  |
| <u>Bulk Solvent Manifold</u>  |                    |  |
|  | <b>01-04200-01</b> | Bulk Solvent Manifold<br>Supplied with 4 meters of 2 mm ID connecting tubing, 10 x Inline Filters<br>Barcode             |



[www.eprep-analytical.com/eprep](http://www.eprep-analytical.com/eprep)

EPREP Pty Ltd  
14/35 Dunlop Rd  
Mulgrave VIC 3170  
AUSTRALIA  
(T) +61-(0)3-9574 3600  
(E) [info@eprep.com.au](mailto:info@eprep.com.au)

**ePREP**<sup>®</sup>  
automation for every laboratory

ePrep® is a registered trademark of ePrep Pty Ltd  
μSPEed® is a registered trademark of ePrep Pty Ltd  
eZy-Connect™ is a trademark of ePrep Pty Ltd  
Other trademarks are trademarks of their respective owner

Publication No.98-15000-03 Feb 2020  
© Copyright 2020, ePrep Pty Ltd