



Analisi di Microplastiche

Analisi per **forma** (LDIR) e per **massa** (Py-GCMS)

Dott. Federico Sacco

Material Market Manager

sacco@srainstruments.com +39 346/3715661


What are Microplastics?

- Microplastics are defined as polymer particles **<5mm**. They are now recognized as an **emerging contaminant** in the environment, however current analytical methods for their **identification** and **quantification** are not fully developed.
- Of most concern are the particles **<50 μ m** as they are the ones that are most likely to be **end up in our food chain**.



Measurement Challenge


- Microplastics are difficult to measure, not only due to their obviously small size, but also because the chosen technique needs to be able to identify a wide variety of polymers existing in a wide variety of different shapes and forms e.g. fibres, fragments and microbeads.
- There is no single “perfect” technique, so a combination of approaches is needed.



Measurement Challenge


Microplastics finds have been reported as number of particles

- Particle counts are not a conserved base quantity – particles are brittle



Different size ranges are used (e.g. 330-5000 μm versus 10-500 μm)

Particle shape is ignored – size is reported by one number only, typically called the “particle size”



All these particles would be treated as having the same size

Particle numbers + sizes are important as impacts increase with decreasing size
BUT
The mass of microplastics is needed to assess and compare sources and contamination levels

Forma



Spettroscopia LDIR

- Dimensioni delle particelle
- Distribuzione morfologica
- Forma (fibre, sferiche, cilindriche,...)
- UDM: n° di particelle/ kg o L di campione


Massa




Pirolisi GCMS

- “Quanta massa di microplastiche ho”
- Analisi quantitative (MP e additivi, std interno o esterno)
- UDM: ug/ kg o L di campione

Forma




Massa



Forma



Massa




Forma



Spettroscopia LDIR

- Dimensioni delle particelle
- Distribuzione morfologica
- Forma (fibre, sferiche, cilindriche,...)
- UDM: n° di particelle/ kg o L di campione





Analytical technique	Shape info	Chemical ID	MP _{Num}	MP _{Mass}
Optical microscopy; Fluorescence microscopy	● ●	● ●	● ●	● ●
ATR-FTIR	N/A	●	N/A	● *
μFT-IR-Imaging	●	●	●	●
NIR, Hy-Spec-Imaging	●	●	●	●
μRaman;	●	●	●	●
Pyr-GC-MS	●	●	●	●
TED-GC-MS	●	●	●	●

* Using a balance to weigh particles


Laser Direct Infrared (LDIR) Spectroscopy

- New instrument architecture
- Bright infrared laser source with high power density (1000x more than a Globar)
- Proprietary Agilent quantum cascade laser (QCL) technology
- Rapidly tunable across the mid-infrared ($1800 - 975\text{cm}^{-1}$)
- Laser spot size at sample is $5.5 - 10\mu\text{m}$ diameter

- Focus all laser power onto a particle
- One second per spectrum
- 0.5cm^{-1} spectral resolution



Coherent laser light




LDIR 8700- QCL Setup


The LDIR is **NOT** an FTIR. It doesn't have an interferometer.

The laser wavelength is selected using a monochromator (same type as we use in the Cary 3500) so it's a form of "dispersive IR" spectrometer.

MCT detector doesn't use LN2 – it's electrically cooled – so always ready to use




LDIR Microplastics Analysis Workflow




Particles	Identifications	Statistics
Diameter	Depth	Quality
51.34 μm	13.00 μm	0.829
1	Polyethylene Terephthalate (PET)	
34.02 μm	6.00 μm	0.878
2	Polypropylene (PP)	
32.11 μm	4.00 μm	0.822
3	Poly(vinyl chloride), carboxylated	
28.12 μm	9.00 μm	0.839
4	Polyamide (PA)	
26.44 μm	6.00 μm	0.876
5	Polyacrylic esters	
25.56 μm	17.00 μm	0.911
6		

Particles	Identifications	Statistics
Polyamide (PA)		43.6%
Polyacrylic esters		10.3%
Polystyrene (PS)		10.3%
Poly(vinyl chloride), carboxylated		7.7%
Polymethyl methacrylate (PMMA)		5.1%
Polyethylene Terephthalate (PET)		5.1%
Cellulose		5.1%
Polypropylene (PP)		2.6%
Polyacrylamide		2.6%
Polycaprolactone		2.6%
Poly(vinyl chloride)		2.6%
Cellulose		2.6%




How do we achieve such great resolution?


Oversampling




10µm Pixel Resolution



5µm Pixel Resolution




3µm Pixel Resolution



1µm Pixel Resolution

- **Polymer mix**


- Blend of nine common microplastics ground using a cryomill and filtered to be between 10-100 μm
 - Polyethylene terephthalate (PET)
 - Polyethylene (PE)
 - Polyvinyl chloride (PVC)
 - Polypropylene (PP)
 - Polystyrene (PS)
 - Polytetrafluoroethylene (PTFE)
 - Polymethyl methacrylate (PMMA)
 - Nylon 6,6 (PA)
 - Styrene butadiene rubber (SBR)
- Particles were suspended in ethanol, sonicated, and aliquoted onto a slide for analysis.







Samples courtesy of
Nikki van Alst / Prof. Jes Vollertsen
Aalborg University, Denmark


The LDIR can handle a wider range of particle sizes – up to 1000 μm (1mm) - in a single analysis, so less prep is needed and faster overall measurement time.

Particle breakdown and statistical analysis



3 mm × 3 mm; 39 particles; 5 minutes





- Both large (1 mm) and small (10 μm) particles can be measured in the same analysis
- Bright laser light is able to transmit through thick particles and ID them

Analisi per massa: Py-GCMS



Massa




Development of analytical method for MPs

**Develop an analytical method and analytical system
for microplastics analysis using Py-GC/MS**




Challenges:

1. Target plastic types
2. Secondary reactions among pyrolyzates (acid, alcohol, amine and unsaturated compounds)
3. Optimization of a separation column for various pyrolyzates
4. Qualification and quantification of plastics (Search software)
5. Matrices difficulties with environmental samples



Schematic diagram of Multi-Shot pyrolyzer EGA/PY-3030D



Target areas

1. Environmental

- MPs in the sea, rivers, lakes, soil, sediments
 - ✓ Waste plastics, Coated fertilizers, Tire wear particles etc.
- MPs in organisms
 - ✓ Planktons, Fish, Waterfowl etc.
- MPS in atmosphere
 - ✓ Tire wear particles etc.



2. Food/drink

- MPs in drinks
 - ✓ Water, Beer, Juice etc.
- MPs in food
 - ✓ Salt, Seafood etc.



3. Others

- Persistent organic pollutants (POPs) on/in MPs
- MPs from packaging materials in food
- MPs from Paint and coating materials in rivers



F-Search MPs 2.0



General product information

F-Search MPs 2.0 is a searching software that assists in the **qualitative** and **quantitative** analysis of microplastics by Pyrolysis(Py-)GC/MS method using Frontier's recommended analytical system and analytical conditions.

We developed an analytical method and configured a recommended Py-GC/MS system for microplastics analysis, so it is easy to start microplastics analysis even for those who are using Py-GC/MS for the first time.



Product features

* Japan patent: #6683335

1) Identify Microplastics polymers by a sophisticated search algorithm*

2) Automated creation of calibration curves and quick quantification

Create calibration curves of the registered polymers (ISTD or ESTD) automatically based on the analytical results from the reference polymer mixture. Then, F-Search MPs 2.0 automatically performs quantification, and the results are instantly reported for the environmental MPs.

3) Cover twelve commonly used polymers


It enables to analyze twelve commonly used polymers (e.g. polyethylene, polypropylene)

4) Enable to create user's own library


It has a feature to establish user's own library depending on the user's interest.

Comparison of database and searching algorithm

F-Search ver. 3.6 Styrene monomer, dimer, and trimer are selected as the characteristic pyrolyzates, and the mass spectrum obtained by Int-Sum for PS is registered in the database library. For qualification, detect peaks with a certain intensity, create Int-Sum mass spectrum for the sample, and search it on the database library.




F-Search MPs 2.0 Styrene monomer is not included as a characteristic pyrolyzate because it can be derived from other styrenic polymer, so the Int-Sum mass spectrum of characteristic ions of dimer and trimer is registered in the database. For qualification, these two peaks and characteristic ions are automatically detected from RI and search it on the database library. A styrene trimer is used for quantification.




Search algorithm of F-Search MPs 2.0

A: Registered polymer in database library (e.g. PS)


Database library




Int-Sum mass spectrum of SS and SSS




B: Unknown sample (Search for PS)




Step 1 :
Detect EIC of characteristic ions of pyrolyzates



Step 2:
Obtain mass spectrum based on Step 1





Step 3:
Sum up mass spectrum of A and B



Step 4:
Calculate matching % of mass spectrum between unknown sample and registered polymer in database library

Analytical procedure





0.21
Run Time: 90.00

EG A/PY-3030D Control - EGA_legno.mtd

File View Tools Help

Direct EGA Analysis

START STOP

Pyrolyzer
 Furnace

Step	Initial (°C)	Initial (min)	Rate (°C/min)	Final (°C)	Final (min)	Total (min)
1st	50	0.00	20	200	0.00	7.50
2nd			8	500	0.00	45.00
3rd			20	700	0.00	55.00
4th			0	40	0.00	

Interface Upper Temp. 300 °C Auto Manual

Accessories
 SS MJT AS CGS UV Settings

Monitor

Not Ready
 Run


Actual Set
 Furnace 50.4 °C 200 °C
 Interface 150.1 °C 300 °C

Elapsed Set
 Time 0.16 min 55.00 min

Elapsed Time

Sample analysis requires a various fases procedure in order to obtain relevant and consistent results

**THIS IS INDIPENDENT FROM THE INSTRUMENT CHOSEN FOR
MPs ANALYSIS**



I vari step di pretrattamento del campione si possono semplificare o rimuovere a seconda della tipologia di matrice che andiamo ad analizzare

Operational process by F-Search MPs 2.0

Step 1 Creating a plot of retention index (RI) and retention time (RT)

RI is created based on the pyrolyzates of Polystyrene using MPs Calibration Standard

Step 2 Creating a calibration file (QFL) (e.g. 3 points, $n=5$)

Calibration file is created based on pyrograms of MPs Calibration Standard


Step 3 Qualitative/quantitative analysis of MPs

Search results are displayed

Step 4 Interpretation of results

Creating a plot of retention index (RI)

- RI is created based on the pyrolyzates of Polystyrene in MP Calibration Standard



Load a pyrogram data of MP Calibration Standard

The 'RI Calc From PS' dialog box contains the following data:

	Retention Index	Retention Time	
		Current	Detected
Toluene:	787	5.884	5.658 [min]
Styrene monomer:	917	7.332	7.359 [min]
Styrene dimer:	1799	13.498	14.508 [min]
Styrene trimer:	2539	17.543	21.344 [min]

Buttons at the bottom: Default RI, OK, Update, Cancel.

Auto-detection of RI

Creating a calibration file (1)

- Select "Add Rep. Files" and upload data files and input sample amount of each polymer.
- Then select Index level (different concentration) "2" and upload data files. You may use multiplier to automatically input sample amount of each polymer.

Max number of files: 10

Cal. Levels: 3 Apply

Index Level: 1

Multiplier: Apply

Polymer	Amount [ug]
PE	160
PP	40
PS	8
ABS	16
SBR	16
PMMA	8
PC	4
PVC	40
PU	10
PET	16
N6	5
N66	18

File name

Ref_CaCO3_1of1_PY_TIAR2101

Ref_CaCO3_1of1_PY_TIAR2101

Ref_CaCO3_1of1_PY_TIAR2101

Ref_CaCO3_1of1_PY_TIAR2101

Ref_CaCO3_1of1_PY_TIAR2101

Polymers with zero set cannot be quantified.

Next

Load data file and input sample amount

Max number of files: 10

Cal. Levels: 3 Apply

Index Level: 2

Multiplier: 0.5 Apply

Polymer	Amount [ug]
PE	80.0
PP	20.0
PS	4.00
ABS	8.00
SBR	8.00
PMMA	4.00
PC	2.00
PVC	20.0
PU	5.00
PET	8.00
N6	2.50
N66	9.00

File name

Ref_CaCO3_1of2_PY...

Ref_CaCO3_1of2_PY...

Ref_CaCO3_1of2_PY...

Ref_CaCO3_1of2_PY...

Ref_CaCO3_1of2_PY...

Polymers with zero set cannot be quantified.

Next

Load data file of different conc.

Max number of files: 10

Cal. Levels: 3 Apply

Index Level: 3

Multiplier: 0.1 Apply

Polymer	Amount [ug]
PE	16.0
PP	4.00
PS	0.80
ABS	1.60
SBR	1.60
PMMA	0.80
PC	0.40
PVC	4.00
PU	1.00
PET	1.60
N6	0.50
N66	1.80

File name

Ref_CaCO3_1of10_PY_TIAR21012505...

Ref_CaCO3_1of10_PY_TIAR21012506...

Ref_CaCO3_1of10_PY_TIAR21012507...

Ref_CaCO3_1of10_PY_TIAR21012508...

Ref_CaCO3_1of10_PY_TIAR21012509...

Polymers with zero set cannot be quantified.

Next

Cancel

Load data file of different conc.

Creating a calibration file (2)

- Automatically creating a calibration file. If there is any problem with peak shape, peak area, peak shape and peak itself, correct it by clicking and dragging a baseline.
- Once correction is completed, select “Save”. A calibration file is created (qfl. file) .

Create Quantification File (2/2)

① Select conc. Index Level: 1 [Delete]

② Select data file # Rep. File: 3/5 [Delete] Ref_CaCO3_1of10_F [<< >>]

Switch a chart of EIC, Peak Area, Peak Height
Chart: EIC Area Height

Poly...	P...	m...	RT [min]	Area	Avr Area	Height	Avr Heig...	Amount ...	Ty...
PE	C...	82	16.349	45185	44479	17566	18200	16.0	BB
PP	C9	1...	6.463	24440	23236	11348	11151	4.00	BB
PS	S...	91	21.332	100438	105693	26153	27320	0.80	BB
ABS	S...	1...	18.114	28660	28818	10410	10455	1.60	BB
SBR	V...	54	6.516	25093	21704	11318	10737	1.60	VIB
PMMA	M...	1...	4.816	85079	86276	32773	32842	0.80	BB
PC	IPP	1...	11.301	114177	115699	63689	63985	0.40	BB
PVC	N...	1...	10.569	74079	74658	41322	41700	4.00	BB
PU	M...	1...	18.008	181658	185529	59009	61178	1.00	BB
PET	BP	1...	14.058	9367	9130	4719	4641	1.60	BB
N6	C...	1...	11.271	31965	31517	16116	15931	0.50	BB
N66	CP	84	6.211	83689	79384	37416	33751	1.80	BB

③ Select a marker peak of each polymer


Enlarged display of EIC. Editable by clicking and dragging a baseline

Save Copy Table

Searching result (1/3) display overview

Results of Qual/Quant*

Mass spectrum and peak shapes



Qual : Qualification
Quant: Quantification

Searching result (2/3) display overview

The screenshot displays the 'MP Search' window with the following components and annotations:

- Threshold of Match quality:** Points to the 'Display Threshold' field, which is set to 0.0 %.
- Library Match Quality (%):** Points to the 'Prob.' column in the results table.
- Plastic Name:** Points to the 'Polymer' column in the results table.
- Quant Result (µg):** Points to the 'Qty.' field, which is set to 0.0 ug.
- Relative Amount (%):** Points to the 'Ratio [%]' column in the results table.
- Threshold of Quant result:** Points to the 'Area' column in the results table.
- Peak area:** Points to the 'Area' column in the results table.
- Retention time:** Points to the 'RT [min]' column in the results table.
- LOQ (calibration, user):** Points to the 'LOQ [...]' column in the results table.
- Bar graph of Relative Amount (%):** Points to the bar chart at the bottom of the window, which shows the relative amounts of various polymers.


Polymer	Prob. [...]	Qty [ug]	Ratio [%]	Area	RT [min]	LOQ [...]
PE	99.9	11.27	49.3	31420	16.36	
PVC	96.6	7.933	34.7	146285	10.57	
SBR	23.9	0.855	3.74	12294	7.07	
PP	92.1	0.724	3.17	4116	6.46	
PS	99.6	0.577	2.53	75144	21.33	
PET	88.6	0.490	2.14	2574	13.92	
PMMA	100.0	0.367	1.61	39050	4.82	
PU	99.9	0.305	1.34	81556	18.01	
ABS	70.5	0.152	0.67	2697	18.02	
N66	88.0	0.111	0.48	4847	6.62	
N6	81.3	0.061	0.27	3745	11.50	
PC	70.8	0.018	0.08	5019	11.24	
(100)						

Qual: Qualification
Quant: Quantification


It is noted that names of plastics in searching results suggest **a polymer group**.

Searching result (3/3) Mass spectrum/Peak shape


Marker peak of Sample




Marker peak of MPs Calibration Standard for calibration



Mass Spectrum of Sample peak




Mass Spectrum of MPs Calibration Standard for calibration



- Each polymer has its own marker peaks and ions

Sample prep for calibration procedure


Put 4 mg of MP Calibration Standard into a sample cup




HR202i (A&D)
< Min: 0.01 mg >



The amount of each polymer in the 4 mg of reference mixture (PE: 160 µg, PP: 40 µg, PS: 8 µg, ABS: 16 µg, SBR: 16 µg, PMMA: 8 µg, PC: 4 µg, PVC: 40 µg, PU: 4 µg, PET: 16 µg, N6: 5 µg, N66: 18 µg, diluent: 3.8 mg)



Summary & Questions ?



Micro Plastic Solution !

- ① Multi-Shot Pyrolyzer
- ② Auto-Shot Sampler
- New** ③ F-Search MPs 2.0
- New** ④ Packed GC Glass Insert
- New** ⑤ UAMP column kit
- ⑥ Vent-free GC/MS adapter
- New** ⑦ Microplastics Calibration Standard
- New** ⑧ Consumable parts set for MPs

UA: Ultra ALLOY®

- Combination with Single Q GC/MS system.
- Interpretation results by F-Search MPs 2.0 are optimized by Py-GC/MS system above.

