


Head Clamp Assembly

For Agilent GC's



FEATURES

- Estimated MTBF > 7 years
- Field Replaceable Module
- Can modulate 2 loops in 1 Column Holder (“buy 1 thermal modulator, get 1 free”)
- Modulates 2 parallel Column Sets
- Hot Jet Gas Fitting: ¼” Swagelok
- Set-and-Forget Installation

SPECIFICATIONS

Max Heater Power (Watts)	80
Max Voltage DC	40
Sensor Type	Platinum Resistor
Sensor Resistance (Ohms)	~ 101
Power Source: Available Aux Thermal Zone	-

Max Heater Temperature (°C)	450
Gas Consumption (L/min)	10 to 15
Typical Gas	Dry N2
Main N2 Supply Pressure (psi)	100
Valve N2 Supply Pressure (psi)	40